

Biotinylated Anti-bovine IL-2 Polyclonal Antibody

KPAB0206

Description

This Biotinylated Anti-bovine IL-2 Polyclonal Antibody has been validated for ELISA, ELISpot applications.

Product Information

SKU:	KPAB0206
Concentration:	1.0 mg/mL
Conjugated:	Biotin
Applications:	ELISA, ELISpot

Antibody Data

Host:	Rabbit
Reactivity:	Bovine
Immunogen:	Recombinant Bovine IL-2
Molecular Weight:	15.5 kDa (calculated) for bovine IL-2
Purification:	Antigen-affinity Purification
Formulation:	1X PBS buffer containing 0.09% NaN ₃
Cross-Reactivity:	Ovine IL-2: Weak
Homology:	Bison bison bison (American buffalo), Bos gaurus (gaur), Bos indicus (zebu), Bos javanicus (banteng), Bos mutus (wild yak), Bos taurus (cattle), Bubalus bubalis (water buffalo)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.