

Biotinylated Anti-feline TNF alpha Polyclonal Antibody

KPAB0209

Description

This Biotinylated Anti-feline TNF alpha Polyclonal Antibody has been validated for ELISA, ELISpot applications.

Product Information

SKU:	KPAB0209
Concentration:	1.0 mg/mL
Conjugated:	Biotin
Applications:	ELISA, ELISpot

Antibody Data

Host:	Rabbit
Reactivity:	Feline
Immunogen:	Recombinant Feline TNF alpha
Molecular Weight:	17.2 kDa (calculated) for feline TNF- α
Purification:	Antigen-affinity Purification
Formulation:	1X PBS buffer containing 0.09% NaN ₃
Cross-Reactivity:	Canine TNF-A: Strong; Caprine TNF-A: Weak; Cynomolgus TNF-A: Strong; Equine TNF-A: Weak; Human TNF-A: Moderate; Swine TNF-A: Weak
Homology:	Felis catus (domestic cat), Prionailurus bengalensis (leopard cat), Prionailurus viverrinus (fishing cat)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.