

Anti-bovine IL-1 beta Polyclonal Antibody

KPAB0210

Description

This Anti-bovine IL-1 beta Polyclonal Antibody has been validated for ELISA, ELISpot, FC, Neutralization, WB applications.

Product Information

| | |
|-----------------------|--|
| SKU: | KPAB0210 |
| Concentration: | 1.0 mg/mL |
| Conjugated: | Unconjugated |
| Applications: | ELISA, ELISpot, FC, Neutralization, WB |

Antibody Data

| | |
|--------------------------|--|
| Host: | Rabbit |
| Reactivity: | Bovine |
| Immunogen: | Recombinant Bovine IL-1 beta |
| Molecular Weight: | 17.7 kDa (calculated) for bovine IL-1 β |
| Purification: | Antigen-affinity Purification |
| Formulation: | 1X PBS buffer containing 0.09% NaN ₃ |
| Cross-Reactivity: | Caprine IL-1B: Strong; Ovine IL-1B: Moderate |
| Homology: | Bos indicus (zebu), Bos javanicus (banteng), Bos taurus (cattle) |

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.