

Biotinylated Anti-chicken IL-6 Polyclonal Antibody

KPAB0213

Description

This Biotinylated Anti-chicken IL-6 Polyclonal Antibody has been validated for ELISA, ELISpot applications.

Product Information

SKU:	KPAB0213
Concentration:	1.0 mg/mL
Conjugated:	Biotin
Applications:	ELISA, ELISpot

Antibody Data

Host:	Rabbit
Reactivity:	Chicken
Immunogen:	Recombinant Chicken IL-6
Molecular Weight:	22.0 kDa (calculated) for chicken IL-6
Purification:	Antigen-affinity Purification
Formulation:	1X PBS buffer containing 0.02% NaN ₃
Cross-Reactivity:	Dolphin IL-6: Weak
Homology:	Gallus gallus (chicken)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.