

## Biotinylated Anti-equine Erythropoietin Polyclonal Antibody

**KPAB0226**

### Description

---

This Biotinylated Anti-equine Erythropoietin Polyclonal Antibody has been validated for ELISA, ELISpot applications.

### Product Information

---

**SKU:** KPAB0226

**Conjugated:** Biotin

**Applications:** ELISA, ELISpot

### Antibody Data

---

**Host:** Rabbit

**Reactivity:** Equine

**Immunogen:** Recombinant Horse Epo

**Purification:** Antigen-affinity Purification

**Homology:** Equus caballus (horse), Equus przewalskii (Przewalski's horse)

## Preparation & Storage

**Storage:** Stable for up to twelve months from date of receipt at 2-8°C.

**Shipping:** Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.