

Anti-dolphin IL-2 Polyclonal Antibody

KPAB0231

Description

This Anti-dolphin IL-2 Polyclonal Antibody has been validated for ELISA, ELISpot, FC, Neutralization, WB applications.

Product Information

| | |
|-----------------------|--|
| SKU: | KPAB0231 |
| Concentration: | 1.0 mg/mL |
| Conjugated: | Unconjugated |
| Applications: | ELISA, ELISpot, FC, Neutralization, WB |

Antibody Data

| | |
|--------------------------|---|
| Host: | Rabbit |
| Reactivity: | Dolphin |
| Immunogen: | Recombinant Dolphin IL-2 |
| Molecular Weight: | 15.5 kDa (calculated) for dolphin IL-2 |
| Purification: | Antigen-affinity Purification |
| Formulation: | 1X PBS buffer containing 0.09% NaN3 |
| Homology: | Globicephala melas (long-finned pilot whale), Orcinus orca (killer whale), Phocoena sinus (vaquita), Tursiops truncatus (common bottlenose dolphin) |

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.