

Anti-equine CCL11 Polyclonal Antibody

KPAB0233

Description

This Anti-equine CCL11 Polyclonal Antibody has been validated for ELISA, ELISpot, FC, Neutralization, WB applications.

Product Information

SKU:	KPAB0233
Concentration:	1.0 mg/mL
Conjugated:	Unconjugated
Applications:	ELISA, ELISpot, FC, Neutralization, WB

Antibody Data

Host:	Rabbit
Reactivity:	Equine
Immunogen:	Recombinant Horse CCL11
Molecular Weight:	9.0 kDa (calculated) for equine CCL11
Purification:	Antigen-affinity Purification
Formulation:	1X PBS buffer containing 0.09% NaN3
Cross-Reactivity:	Bovine CCL11: Moderate; Caprine CCL11: Moderate; Guinea Pig CCL11: Weak
Homology:	Equus asinus (ass), Equus caballus (horse), Equus przewalskii (Przewalski's horse), Equus quagga (plains zebra)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.