

## Biotinylated Anti-equine CCL11 Polyclonal Antibody

**KPAB0234**

### Description

---

This Biotinylated Anti-equine CCL11 Polyclonal Antibody has been validated for ELISA, ELISpot applications.

### Product Information

---

**SKU:** KPAB0234  
**Conjugated:** Biotin  
**Applications:** ELISA, ELISpot

### Antibody Data

---

**Host:** Rabbit  
**Reactivity:** Equine  
**Immunogen:** Recombinant Horse CCL11  
**Purification:** Antigen-affinity Purification  
**Homology:** Equus asinus (ass), Equus caballus (horse), Equus przewalskii (Przewalski's horse), Equus quagga (plains zebra)

## Preparation & Storage

**Storage:** Stable for up to twelve months from date of receipt at 2-8°C.

**Shipping:** Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.