

Biotinylated Anti-bovine IL-8 Polyclonal Antibody

KPAB0235

Description

This Biotinylated Anti-bovine IL-8 Polyclonal Antibody has been validated for ELISA, ELISpot applications.

Product Information

SKU:	KPAB0235
Concentration:	1.0 mg/mL
Conjugated:	Biotin
Applications:	ELISA, ELISpot

Antibody Data

Host:	Rabbit
Reactivity:	Bovine
Immunogen:	Recombinant Bovine IL-8
Molecular Weight:	9.1 kDa (calculated) for bovine IL-8
Purification:	Antigen-affinity Purification
Formulation:	1X PBS buffer containing 0.09% NaN ₃
Cross-Reactivity:	Canine IL-8: Moderate; Dolphin IL-8: Moderate; Equine IL-8: Weak; Ovine IL-8: Moderate; Swine IL-8: Weak
Homology:	Bos indicus (zebu), Bos taurus (cattle)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.