

Biotinylated Anti-bovine IL-8 Polyclonal Antibody

KPAB0237

Description

This Biotinylated Anti-bovine IL-8 Polyclonal Antibody has been validated for ELISA, ELISpot applications.

Product Information

SKU: KPAB0237
Conjugated: Biotin
Applications: ELISA, ELISpot

Antibody Data

Host: Goat
Reactivity: Bovine
Immunogen: Recombinant Bovine IL-8
Purification: Antigen-affinity Purification
Homology: Bos indicus (zebu), Bos taurus (cattle)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.