

Anti-cynomolgus monkey IL-6 Polyclonal Antibody

KPAB0238

Description

This Anti-cynomolgus monkey IL-6 Polyclonal Antibody has been validated for ELISA, ELISpot, FC, Neutralization, WB applications.

Product Information

SKU:	KPAB0238
Concentration:	1.0 mg/mL
Conjugated:	Unconjugated
Applications:	ELISA, ELISpot, FC, Neutralization, WB

Antibody Data

Host:	Rabbit
Reactivity:	Cynomolgus Monkey
Immunogen:	Recombinant Cynomolgus Monkey IL-6
Molecular Weight:	21.0 kDa (calculated) for cynomolgus monkey IL-6
Purification:	Antigen-affinity Purification
Formulation:	1X PBS buffer containing 0.09% NaN ₃
Cross-Reactivity:	Human IL-6: Strong
Homology:	Macaca fascicularis (crab-eating macaque), Macaca mulatta (Rhesus monkey), Macaca nemestrina (pig-tailed macaque)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.