

## Anti-swine IL-1 alpha Polyclonal Antibody

### KPAB0240

#### Description

---

This Anti-swine IL-1 alpha Polyclonal Antibody has been validated for ELISA, ELISpot, FC, Neutralization, WB applications.

#### Product Information

---

<b>SKU:</b>	KPAB0240
<b>Concentration:</b>	1.0 mg/mL
<b>Conjugated:</b>	Unconjugated
<b>Applications:</b>	ELISA, ELISpot, FC, Neutralization, WB

#### Antibody Data

---

<b>Host:</b>	Rabbit
<b>Reactivity:</b>	Swine
<b>Immunogen:</b>	Recombinant Porcine IL-1 alpha
<b>Molecular Weight:</b>	18.1 kDa (calculated) for swine IL-1 $\alpha$
<b>Purification:</b>	Antigen-affinity Purification
<b>Formulation:</b>	1X PBS buffer containing 0.09% NaN <sub>3</sub>
<b>Cross-Reactivity:</b>	Bovine IL-1A: Moderate; Caprine IL-1A: Moderate; Equine IL-1A: Moderate; Feline IL-1A: Weak
<b>Homology:</b>	Phacochoerus africanus (Common warthog), Sus scrofa (pig)

## Preparation & Storage

**Storage:** Stable for up to twelve months from date of receipt at 2-8°C.

**Shipping:** Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.