

Biotinylated Anti-human IFN gamma Polyclonal Antibody

KPAB0243

Description

This Biotinylated Anti-human IFN gamma Polyclonal Antibody has been validated for ELISA, ELISpot applications.

Product Information

SKU:	KPAB0243
Concentration:	1.0 mg/mL
Conjugated:	Biotin
Applications:	ELISA, ELISpot

Antibody Data

Host:	Goat
Reactivity:	Human
Immunogen:	Recombinant Human IFN gamma
Molecular Weight:	16.8 kDa (calculated) for human IFN- η
Purification:	Antigen-affinity Purification
Formulation:	1X PBS buffer containing 0.09% NaN ₃
Cross-Reactivity:	Cynomolgus IFN-H: Strong
Homology:	Gorilla gorilla (western gorilla), Homo sapiens (human)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.