

## Anti-bovine CD40 Ligand Polyclonal Antibody

**KPAB0244**

### Description

---

This Anti-bovine CD40 Ligand Polyclonal Antibody has been validated for ELISA, ELISpot, FC, Neutralization, WB applications.

### Product Information

---

<b>SKU:</b>	KPAB0244
<b>Concentration:</b>	1.0 mg/mL
<b>Conjugated:</b>	Unconjugated
<b>Applications:</b>	ELISA, ELISpot, FC, Neutralization, WB

### Antibody Data

---

<b>Host:</b>	Rabbit
<b>Reactivity:</b>	Bovine
<b>Immunogen:</b>	Recombinant Bovine CD40 Ligand
<b>Molecular Weight:</b>	16.1 kDa (calculated) for bovine CD40LG
<b>Purification:</b>	Antigen-affinity Purification
<b>Formulation:</b>	1X PBS buffer containing 0.09% NaN <sub>3</sub>
<b>Cross-Reactivity:</b>	Canine CCL4: Moderate; Equine CCL4: Strong; Feline CCL4: Strong; Human CCL4: Moderate; Swine CCL4: Strong
<b>Homology:</b>	Bison bison bison (American buffalo), Bos mutus (wild yak), Bos taurus (cattle)

## Preparation & Storage

**Storage:** Stable for up to twelve months from date of receipt at 2-8°C.

**Shipping:** Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.