

Biotinylated Anti-bovine CD40 Ligand Polyclonal Antibody

KPAB0245

Description

This Biotinylated Anti-bovine CD40 Ligand Polyclonal Antibody has been validated for ELISA, ELISpot applications.

Product Information

SKU:	KPAB0245
Concentration:	1.0 mg/mL
Conjugated:	Biotin
Applications:	ELISA, ELISpot

Antibody Data

Host:	Rabbit
Reactivity:	Bovine
Immunogen:	Recombinant Bovine CD40 Ligand
Molecular Weight:	16.1 kDa (calculated) for bovine CD40LG
Purification:	Antigen-affinity Purification
Formulation:	1X PBS buffer containing 0.09% NaN ₃
Cross-Reactivity:	Canine CCL4: Moderate; Equine CCL4: Strong; Feline CCL4: Strong; Human CCL4: Moderate; Swine CCL4: Strong
Homology:	Bison bison bison (American buffalo), Bos mutus (wild yak), Bos taurus (cattle)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.