

## Anti-equine IL-12 p40 Polyclonal Antibody

### KPAB0250

#### Description

---

This Anti-equine IL-12 p40 Polyclonal Antibody has been validated for ELISA, ELISpot, FC, Neutralization, WB applications.

#### Product Information

---

<b>SKU:</b>	KPAB0250
<b>Concentration:</b>	1.0 mg/mL
<b>Conjugated:</b>	Unconjugated
<b>Applications:</b>	ELISA, ELISpot, FC, Neutralization, WB

#### Antibody Data

---

<b>Host:</b>	Rabbit
<b>Reactivity:</b>	Equine
<b>Immunogen:</b>	Recombinant Horse IL-12 p40
<b>Molecular Weight:</b>	34.9 kDa (calculated) for equine IL-12 p40
<b>Purification:</b>	Antigen-affinity Purification
<b>Formulation:</b>	1X PBS buffer containing 0.09% NaN3
<b>Cross-Reactivity:</b>	Bovine IL-12 P40: Weak; Chicken IL-12 P40: Weak; Feline IL-12 P40: Moderate; Ovine IL-12 P40: Moderate; Swine IL-12 P40: Moderate
<b>Homology:</b>	Equus caballus (horse)

## Preparation & Storage

**Storage:** Stable for up to twelve months from date of receipt at 2-8°C.

**Shipping:** Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.