

Biotinylated Anti-equine IL-12 p40 Polyclonal Antibody

KPAB0251

Description

This Biotinylated Anti-equine IL-12 p40 Polyclonal Antibody has been validated for ELISA, ELISpot applications.

Product Information

SKU: KPAB0251
Conjugated: Biotin
Applications: ELISA, ELISpot

Antibody Data

Host: Rabbit
Reactivity: Equine
Immunogen: Recombinant Horse IL-12 p40
Purification: Antigen-affinity Purification
Homology: Equus caballus (horse)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.