

Biotinylated Anti-zebrafish IFN gamma 1-1 Polyclonal Antibody

KPAB0255

Description

This Biotinylated Anti-zebrafish IFN gamma 1-1 Polyclonal Antibody has been validated for ELISA, ELISpot applications.

Product Information

SKU: KPAB0255

Conjugated: Biotin

Applications: ELISA, ELISpot

Antibody Data

Host: Rabbit

Reactivity: Zebrafish

Immunogen: Recombinant Zebrafish IFN gamma 1-1

Purification: Antigen-affinity Purification

Homology: Danio rerio (zebrafish)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.