

Anti-zebrafish IL-17 a/f1 Polyclonal Antibody

KPAB0258

Description

This Anti-zebrafish IL-17 a/f1 Polyclonal Antibody has been validated for ELISA, ELISpot, FC, Neutralization, WB applications.

Product Information

SKU:	KPAB0258
Concentration:	1.0 mg/mL
Conjugated:	Unconjugated
Applications:	ELISA, ELISpot, FC, Neutralization, WB

Antibody Data

Host:	Rabbit
Reactivity:	Zebrafish
Immunogen:	Recombinant Zebrafish IL-17 a/f1
Molecular Weight:	14.4 kDa (calculated) for zebrafish IL-17a/f1
Purification:	Antigen-affinity Purification
Formulation:	1X PBS buffer containing 0.09% NaN ₃
Homology:	Danio rerio (zebrafish)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.