

## Biotinylated Anti-zebrafish IL-17 a/f1 Polyclonal Antibody

**KPAB0259**

### Description

---

This Biotinylated Anti-zebrafish IL-17 a/f1 Polyclonal Antibody has been validated for ELISA, ELISpot applications.

### Product Information

---

**SKU:** KPAB0259  
**Conjugated:** Biotin  
**Applications:** ELISA, ELISpot

### Antibody Data

---

**Host:** Rabbit  
**Reactivity:** Zebrafish  
**Immunogen:** Recombinant Zebrafish IL-17 a/f1  
**Purification:** Antigen-affinity Purification  
**Homology:** Danio rerio (zebrafish)

## Preparation & Storage

**Storage:** Stable for up to twelve months from date of receipt at 2-8°C.

**Shipping:** Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.