

## Anti-cynomolgus monkey IL-17A Polyclonal Antibody

**KPAB0260**

### Description

---

This Anti-cynomolgus monkey IL-17A Polyclonal Antibody has been validated for ELISA, ELISpot, FC, Neutralization, WB applications.

### Product Information

---

<b>SKU:</b>	KPAB0260
<b>Concentration:</b>	1.0 mg/mL
<b>Conjugated:</b>	Unconjugated
<b>Applications:</b>	ELISA, ELISpot, FC, Neutralization, WB

### Antibody Data

---

<b>Host:</b>	Rabbit
<b>Reactivity:</b>	Cynomolgus Monkey
<b>Immunogen:</b>	Recombinant Cynomolgus Monkey IL-17A
<b>Molecular Weight:</b>	15.1 kDa (calculated) for cynomolgus monkey IL-17A
<b>Purification:</b>	Antigen-affinity Purification
<b>Formulation:</b>	1X PBS buffer containing 0.09% NaN <sub>3</sub>
<b>Cross-Reactivity:</b>	Canine IL-17A: Strong; Dolphin IL-17A: Weak; Equine IL-17A: Strong; Ferret IL-17A: Strong; Human IL-17A: Strong; Ovine IL-17A: Weak; Rabbit IL-17A: Weak; Rat IL-17A: Weak; Swine IL-17A: Weak
<b>Homology:</b>	Macaca fascicularis (crab-eating macaque), Macaca nemestrina (pig-tailed macaque)

## Preparation & Storage

**Storage:** Stable for up to twelve months from date of receipt at 2-8°C.

**Shipping:** Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.