

Biotinylated Anti-cynomolgus monkey IL-17A Polyclonal Antibody

KPAB0261

Description

This Biotinylated Anti-cynomolgus monkey IL-17A Polyclonal Antibody has been validated for ELISA, ELISpot applications.

Product Information

SKU:	KPAB0261
Concentration:	1.0 mg/mL
Conjugated:	Biotin
Applications:	ELISA, ELISpot

Antibody Data

Host:	Rabbit
Reactivity:	Cynomolgus Monkey
Immunogen:	Recombinant Cynomolgus Monkey IL-17A
Molecular Weight:	15.1 kDa (calculated) for cynomolgus monkey IL-17A
Purification:	Antigen-affinity Purification
Formulation:	1X PBS buffer containing 0.09% NaN ₃
Cross-Reactivity:	Canine IL-17A: Strong; Dolphin IL-17A: Weak; Equine IL-17A: Strong; Ferret IL-17A: Strong; Human IL-17A: Strong; Ovine IL-17A: Weak; Rabbit IL-17A: Weak; Rat IL-17A: Weak; Swine IL-17A: Weak
Homology:	Macaca fascicularis (crab-eating macaque), Macaca nemestrina (pig-tailed macaque)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.