

Biotinylated Anti-bovine IL-10 Polyclonal Antibody

KPAB0263

Description

This Biotinylated Anti-bovine IL-10 Polyclonal Antibody has been validated for ELISA, ELISpot applications.

Product Information

SKU:	KPAB0263
Concentration:	1.0 mg/mL
Conjugated:	Biotin
Applications:	ELISA, ELISpot

Antibody Data

Host:	Rabbit
Reactivity:	Bovine
Immunogen:	Recombinant Bovine IL-10
Molecular Weight:	18.6 kDa (calculated) for bovine IL-10
Purification:	Antigen-affinity Purification
Formulation:	1X PBS buffer containing 0.09% NaN ₃
Cross-Reactivity:	Dolphin IL-10: Strong; Ovine IL-10: Weak; Swine IL-10: Moderate
Homology:	Bos taurus (cattle)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.