

## Anti-cynomolgus monkey CCL11 Polyclonal Antibody

**KPAB0266**

### Description

---

This Anti-cynomolgus monkey CCL11 Polyclonal Antibody has been validated for ELISA, ELISpot, FC, Neutralization, WB applications.

### Product Information

---

<b>SKU:</b>	KPAB0266
<b>Conjugated:</b>	Unconjugated
<b>Applications:</b>	ELISA, ELISpot, FC, Neutralization, WB

### Antibody Data

---

<b>Host:</b>	Rabbit
<b>Reactivity:</b>	Cynomolgus Monkey
<b>Immunogen:</b>	Recombinant Cynomolgus Monkey CCL11
<b>Purification:</b>	Antigen-affinity Purification
<b>Homology:</b>	Macaca fascicularis (crab-eating macaque), Macaca nemestrina (pig-tailed macaque), Mandrillus leucophaeus (drill)

## Preparation & Storage

**Storage:** Stable for up to twelve months from date of receipt at 2-8°C.

**Shipping:** Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.