

Anti-cynomolgus monkey CXCL10 Polyclonal Antibody

KPAB0268

Description

This Anti-cynomolgus monkey CXCL10 Polyclonal Antibody has been validated for WB, ELISA applications.

Product Information

SKU:	KPAB0268
Conjugated:	Unconjugated
Applications:	WB, ELISA

Antibody Data

Host:	Rabbit
Reactivity:	Cynomolgus Monkey
Immunogen:	Recombinant Cynomolgus Monkey CXCL10
Purification:	Antigen-affinity Purification
Homology:	Cercocebus atys (sooty mangabey), Macaca fascicularis (crab-eating macaque), Macaca mulatta (Rhesus monkey), Mandrillus leucophaeus (drill), Papio anubis (olive baboon), Rhinopithecus bieti (black snub-nosed monkey), Rhinopithecus roxellana (golden snub-nosed monkey), Theropithecus gelada (gelada), Trachypithecus francoisi (Francois's langur)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.