

Biotinylated Anti-cynomolgus monkey CXCL10 Polyclonal Antibody

KPAB0269

Description

This Biotinylated Anti-cynomolgus monkey CXCL10 Polyclonal Antibody has been validated for ELISA, ELISpot applications.

Product Information

SKU:	KPAB0269
Concentration:	1.0 mg/mL
Conjugated:	Biotin
Applications:	ELISA, ELISpot

Antibody Data

Host:	Rabbit
Reactivity:	Cynomolgus Monkey
Immunogen:	Recombinant Cynomolgus Monkey CXCL10
Molecular Weight:	8.7 kDa (calculated) for cynomolgus CXCL10
Purification:	Antigen-affinity Purification
Formulation:	1X PBS buffer containing 0.09% NaN ₃
Homology:	Cercocebus atys (sooty mangabey), Macaca fascicularis (crab-eating macaque), Macaca mulatta (Rhesus monkey), Mandrillus leucophaeus (drill), Papio anubis (olive baboon), Rhinopithecus bieti (black snub-nosed monkey), Rhinopithecus roxellana (golden snub-nosed monkey), Theropithecus gelada (gelada), Trachypithecus francoisi (Francois's langur)
Cross-Reactivity:	Human CXCL10: Strong

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.