

## Biotinylated Anti-canine Fas Ligand Polyclonal Antibody

**KPAB0275**

### Description

---

This Biotinylated Anti-canine Fas Ligand Polyclonal Antibody has been validated for ELISA, ELISpot applications.

### Product Information

---

**SKU:** KPAB0275

**Conjugated:** Biotin

**Applications:** ELISA, ELISpot

### Antibody Data

---

**Host:** Rabbit

**Reactivity:** Canine

**Immunogen:** Recombinant Canine Fas Ligand

**Purification:** Antigen-affinity Purification

**Homology:** Canis lupus dingo (dingo), Canis lupus familiaris (dog), Vulpes lagopus (Arctic fox), Vulpes vulpes (red fox)

## Preparation & Storage

**Storage:** Stable for up to twelve months from date of receipt at 2-8°C.

**Shipping:** Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.