

Anti-dolphin IL-13 Polyclonal Antibody

KPAB0278

Description

This Anti-dolphin IL-13 Polyclonal Antibody has been validated for ELISA, ELISpot, FC, Neutralization, WB applications.

Product Information

SKU:	KPAB0278
Concentration:	1.0 mg/mL
Conjugated:	Unconjugated
Applications:	ELISA, ELISpot, FC, Neutralization, WB

Antibody Data

Host:	Rabbit
Reactivity:	Dolphin
Immunogen:	Recombinant Dolphin IL-13
Molecular Weight:	12.4 kDa (calculated) for dolphin IL-13
Purification:	Antigen-affinity Purification
Formulation:	1X PBS buffer containing 0.09% NaN ₃
Homology:	Globicephala melas (long-finned pilot whale), Neophocaena asiaeorientalis asiaeorientalis (Yangtze finless porpoise), Orcinus orca (killer whale), Phocoena sinus (vaquita), Tursiops truncatus (common bottlenose dolphin)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.