

Biotinylated Anti-dolphin IL-13 Polyclonal Antibody

KPAB0279

Description

This Biotinylated Anti-dolphin IL-13 Polyclonal Antibody has been validated for ELISA, ELISpot applications.

Product Information

SKU: KPAB0279
Conjugated: Biotin
Applications: ELISA, ELISpot

Antibody Data

Host: Rabbit
Reactivity: Dolphin
Immunogen: Recombinant Dolphin IL-13
Purification: Antigen-affinity Purification
Homology: Globicephala melas (long-finned pilot whale), Neophocaena asiaeorientalis asiaeorientalis (Yangtze finless porpoise), Orcinus orca (killer whale), Phocoena sinus (vaquita), Tursiops truncatus (common bottlenose dolphin)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.