

## Biotinylated Anti-swine IL-28B Polyclonal Antibody

**KPAB0281**

### Description

---

This Biotinylated Anti-swine IL-28B Polyclonal Antibody has been validated for ELISA, ELISpot applications.

### Product Information

---

**SKU:** KPAB0281  
**Conjugated:** Biotin  
**Applications:** ELISA, ELISpot

### Antibody Data

---

**Host:** Rabbit  
**Reactivity:** Swine  
**Immunogen:** Recombinant Porcine IL-28B  
**Purification:** Antigen-affinity Purification  
**Homology:** Phacochoerus africanus (Common warthog), Sus scrofa (pig)

## Preparation & Storage

**Storage:** Stable for up to twelve months from date of receipt at 2-8°C.

**Shipping:** Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.