

Biotinylated Anti-rabbit IFN gamma Polyclonal Antibody

KPAB0285

Description

This Biotinylated Anti-rabbit IFN gamma Polyclonal Antibody has been validated for ELISA, ELISpot applications.

Product Information

| | |
|-----------------------|----------------|
| SKU: | KPAB0285 |
| Concentration: | 1.0 mg/mL |
| Conjugated: | Biotin |
| Applications: | ELISA, ELISpot |

Antibody Data

| | |
|--------------------------|---|
| Host: | Chicken |
| Reactivity: | Rabbit |
| Immunogen: | Recombinant Rabbit IFN gamma |
| Molecular Weight: | 16.9 kDa (calculated) for rabbit IFN- η |
| Purification: | Antigen-affinity Purification |
| Formulation: | 1X PBS buffer containing 0.09% NaN ₃ |
| Cross-Reactivity: | Human IFN-H: Weak |
| Homology: | Oryctolagus cuniculus (rabbit) |

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.