

Anti-IGF-1 Polyclonal Antibody

KPAB0290

Description

This Anti-IGF-1 Polyclonal Antibody has been validated for ELISA, ELISpot, FC, Neutralization, WB applications.

Product Information

SKU:	KPAB0290
Concentration:	1.0 mg/mL
Conjugated:	Unconjugated
Applications:	ELISA, ELISpot, FC, Neutralization, WB

Antibody Data

Host:	Goat
Reactivity:	Multi-species
Immunogen:	Recombinant IGF-1
Molecular Weight:	7.6 kDa (calculated) for IGF-1
Purification:	Antigen-affinity Purification
Formulation	1X PBS buffer containing 0.09% NaN ₃
Homology:	Ailuropoda melanoleuca (giant panda), Bison bison bison (American Buffalo), Bos indicus (zebu), Bos mutus (wild yak), Bos taurus (cattle), Bubalus bubalis (water buffalo), Camelus bactrianus (Bactrian camel), Camelus ferus (Wild Bactrian camel), Canis lupus familiaris (dog), Capra hircus (goat), Carlito syrichta (Philippine tarsier), Cavia porcellus (domestic guinea pig), Ceratotherium simum simum (southern white rhinoceros), Cercocebus atys (sooty mangabey), Chinchilla lanigera (long-tailed chinchilla), Chlorocebus sabaeus (green monkey), Dipodomys ordii (Ord's kangaroo rat), Equus asinus (ass), Equus caballus (horse), Equus przewalskii (Przewalski's horse), Felis catus (domestic cat), Gorilla gorilla gorilla (western lowland gorilla), Homo sapiens (human), Ictidomys

tridecemlineatus (thirteen-lined ground squirrel), *Leptonychotes weddellii* (Weddell seal), *Lipotes vexillifer* (Yangtze River dolphin), *Macaca fascicularis* (crab-eating macaque), *Macaca mulatta* (Rhesus monkey), *Macaca nemestrina* (pig-tailed macaque), *Mandrillus leucophaeus* (drill), *Manis javanica* (Malayan pangolin), *Marmota marmota marmota* (Alpine marmot), *Microcebus murinus* (gray mouse lemur), *Mustela putorius furo* (domestic ferret), *Nannospalax galili* (Upper Galilee mountains blind mole rat), *Neogale vison* (American mink), *Nomascus leucogenys* (northern white-cheeked gibbon), *Odobenus rosmarus divergens* (Pacific walrus), *Orcinus orca* (killer whale), *Otolemur garnettii* (small-eared galago), *Pan paniscus* (pygmy chimpanzee), *Pan troglodytes* (chimpanzee), *Panthera pardus* (leopard), *Panholops hodgsonii* (chiru), *Papio anubis* (olive baboon), *Physeter catodon* (sperm whale), *Pongo abelii* (Sumatran orangutan), *Pteropus giganteus* (Indian flying fox), *Pteropus vampyrus* (large flying fox), *Rhinolophus sinicus* (Chinese rufous horseshoe bat), *Rhinopithecus bieti* (black snub-nosed monkey), *Rhinopithecus roxellana* (golden snub-nosed monkey), *Sus scrofa* (pig), *Tursiops truncatus* (bottlenose dolphin), *Vicugna pacos* (alpaca)

Preparation & Storage

- Storage:** Stable for up to twelve months from date of receipt at 2-8°C.
- Shipping:** Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.