

Biotinylated Anti-canine CCL3 Polyclonal Antibody

KPAB0301

Description

This Biotinylated Anti-canine CCL3 Polyclonal Antibody has been validated for ELISA, ELISpot applications.

Product Information

SKU: KPAB0301
Conjugated: Biotin
Applications: ELISA, ELISpot

Antibody Data

Host: Rabbit
Reactivity: Canine
Immunogen: Recombinant Canine CCL3
Purification: Antigen-affinity Purification
Homology: Canis lupus familiaris (dog)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.