

Anti-rabbit VEGF-A Polyclonal Antibody

KPAB0308

Description

This Anti-rabbit VEGF-A Polyclonal Antibody has been validated for ELISA, ELISpot, FC, Neutralization, WB applications.

Product Information

SKU:	KPAB0308
Concentration:	1.0 mg/mL
Conjugated:	Unconjugated
Applications:	ELISA, ELISpot, FC, Neutralization, WB

Antibody Data

Host:	Chicken
Reactivity:	Rabbit
Immunogen:	Recombinant Rabbit VEGF-A 121
Molecular Weight:	14.2 kDa (calculated) for Rabbit VEGF-A121
Purification:	Antigen-affinity Purification
Formulation:	1X PBS buffer containing 0.09% NaN ₃
Homology:	Oryctolagus cuniculus (rabbit)
Cross-Reactivity:	Bovine VEGF-A: Moderate; Canine VEGF-A: Strong; Caprine VEGF-A: Moderate; Equine VEGF-A: Moderate; Feline VEGF-A: Strong; Human/Monkey VEGF-A: Strong; Mouse VEGF-A: Moderate; Ovine VEGF-A: Moderate; Rabbit VEGF-A165: Strong; Rat VEGF-A: Moderate; Swine VEGF-A: Strong

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.