

## Biotinylated Anti-rabbit IL-6 Polyclonal Antibody

### KPAB0315

#### Description

---

This Biotinylated Anti-rabbit IL-6 Polyclonal Antibody has been validated for ELISA, ELISpot applications.

#### Product Information

---

|                       |                |
|-----------------------|----------------|
| <b>SKU:</b>           | KPAB0315       |
| <b>Concentration:</b> | 1.0 mg/mL      |
| <b>Conjugated:</b>    | Biotin         |
| <b>Applications:</b>  | ELISA, ELISpot |

#### Antibody Data

---

|                          |   |
|--------------------------|---|
| <b>Host:</b>             | Chicken   |
| <b>Reactivity:</b>       | Rabbit  |
| <b>Immunogen:</b>        | Recombinant Rabbit IL-6                         |
| <b>Molecular Weight:</b> | 24.5 kDa (calculated) for Rabbit IL-6           |
| <b>Purification:</b>     | Antigen-affinity Purification                   |
| <b>Formulation:</b>      | 1X PBS buffer containing 0.09% NaN <sub>3</sub> |
| <b>Homology:</b>         | Oryctolagus cuniculus (rabbit)                  |

## Preparation & Storage

**Storage:** Stable for up to twelve months from date of receipt at 2-8°C.

**Shipping:** Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.