

## Anti-p38 MAPK Polyclonal Antibody

KPAB0317

### Description

---

This Anti-p38 MAPK Polyclonal Antibody has been validated for WB applications.

### Product Information

---

**SKU:** KPAB0317

**Conjugated:** Unconjugated

**Applications:** WB

### Antibody Data

---

**Host:** Rabbit

**Reactivity:** Multi-species

**Immunogen:** AYFAQYHDPDDEPVADPYDC

**Purification:** Antigen-affinity Purification

## Preparation & Storage

**Storage:** Stable for up to twelve months from date of receipt at 2-8°C.

**Shipping:** Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.