

## Anti-PTEN Polyclonal Antibody

**KPAB0318**

### Description

---

This Anti-PTEN Polyclonal Antibody has been validated for WB applications.

### Product Information

---

**SKU:** KPAB0318  
**Conjugated:** Unconjugated  
**Applications:** WB

### Antibody Data

---

**Host:** Rabbit  
**Reactivity:** Multi-species  
**Immunogen:** STSVTPDVSDNEPDHYRYSC  
**Purification:** Antigen-affinity Purification

## Preparation & Storage

**Storage:** Stable for up to twelve months from date of receipt at 2-8°C.

**Shipping:** Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.