

Anti-ovine IL-22 Polyclonal Antibody

KPAB0323

Description

This Anti-ovine IL-22 Polyclonal Antibody has been validated for ELISA, ELISpot, FC, Neutralization, WB applications.

Product Information

SKU:	KPAB0323
Concentration:	0.5 mg/mL
Conjugated:	Unconjugated
Applications:	ELISA, ELISpot, FC, Neutralization, WB

Antibody Data

Host:	Rabbit
Reactivity:	Ovine
Immunogen:	Recombinant Sheep IL-22
Molecular Weight:	18.3 kDa (calculated) for ovine IL-22
Purification:	Antigen-affinity Purification
Formulation:	1X PBS buffer containing 0.09% NaN3
Cross-Reactivity:	Caprine IL-22: Moderate
Homology:	Ovis aries (sheep)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.