

Anti-cynomolgus monkey CCL2 Polyclonal Antibody

KPAB0337

Description

This Anti-cynomolgus monkey CCL2 Polyclonal Antibody has been validated for ELISA, ELISpot, FC, Neutralization, WB applications.

Product Information

SKU:	KPAB0337
Concentration:	1.0 mg/mL
Conjugated:	Unconjugated
Applications:	ELISA, ELISpot, FC, Neutralization, WB

Antibody Data

Host:	Rabbit
Reactivity:	Cynomolgus Monkey
Immunogen:	Recombinant Cynomolgus Monkey CCL2
Molecular Weight:	8.7 kDa (calculated) for cynomolgus monkey CCL2
Purification:	Antigen-affinity Purification
Formulation:	1X PBS buffer containing 0.09% NaN ₃
Homology:	Cercocebus atys (sooty mangabey), Chlorocebus sabaeus (green monkey), Macaca fascicularis (crab-eating macaque), Macaca mulatta (Rhesus monkey), Macaca nemestrina (pig-tailed macaque), Mandrillus leucophaeus (drill), Papio anubis (olive baboon), Pongo pygmaeus (Bornean orangutan), Rhinopithecus bieti (black snub-nosed monkey), Rhinopithecus roxellana (golden snub-nosed monkey)
Cross-Reactivity:	Bovine CCL2: Moderate; Canine CCL2: Strong; Equine CCL2: Strong; Swine CCL2: Strong

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.