

Biotinylated Anti-cynomolgus monkey CCL2 Polyclonal Antibody

KPAB0338

Description

This Biotinylated Anti-cynomolgus monkey CCL2 Polyclonal Antibody has been validated for ELISA, ELISpot applications.

Product Information

SKU: KPAB0338

Conjugated: Biotin

Applications: ELISA, ELISpot

Antibody Data

Host: Rabbit

Reactivity: Cynomolgus Monkey

Immunogen: Recombinant Cynomolgus Monkey CCL2

Purification: Antigen-affinity Purification

Homology: Cercopithecus atys (sooty mangabey), Chlorocebus sabaeus (green monkey), Macaca fascicularis (crab-eating macaque), Macaca mulatta (Rhesus monkey), Macaca nemestrina (pig-tailed macaque), Mandrillus leucophaeus (drill), Papio anubis (olive baboon), Pongo pygmaeus (Bornean orangutan), Rhinopithecus bieti (black snub-nosed monkey), Rhinopithecus roxellana (golden snub-nosed monkey)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.