

Biotinylated Anti-guinea pig IL-17A Polyclonal Antibody

KPAB0340

Description

This Biotinylated Anti-guinea pig IL-17A Polyclonal Antibody has been validated for ELISA, ELISpot applications.

Product Information

SKU: KPAB0340
Conjugated: Biotin
Applications: ELISA, ELISpot

Antibody Data

Host: Rabbit
Reactivity: Guinea Pig
Immunogen: Recombinant Guinea Pig IL-17A
Purification: Antigen-affinity Purification
Homology: Cavia porcellus (domestic guinea pig)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.