

Biotinylated Anti-ferret IFN gamma Polyclonal Antibody

KPAB0342

Description

This Biotinylated Anti-ferret IFN gamma Polyclonal Antibody has been validated for ELISA, ELISpot applications.

Product Information

SKU: KPAB0342
Conjugated: Biotin
Applications: ELISA, ELISpot

Antibody Data

Host: Rabbit
Reactivity: Ferret
Immunogen: Recombinant Ferret IFN gamma
Purification: Antigen-affinity Purification
Homology: Mustela putorius furo (domestic ferret)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.