

Biotinylated Anti-human EGF Polyclonal Antibody

KPAB0344

Description

This Biotinylated Anti-human EGF Polyclonal Antibody has been validated for ELISA, ELISpot applications.

Product Information

SKU: KPAB0344
Conjugated: Biotin
Applications: ELISA, ELISpot

Antibody Data

Host: Rabbit
Reactivity: Human
Immunogen: Recombinant Human EGF
Purification: Antigen-affinity Purification
Homology: Gorilla gorilla gorilla (western lowland gorilla), Homo sapiens (human), Pan paniscus (pygmy chimpanzee), Pan troglodytes (chimpanzee)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.