

Biotinylated Anti-human IL-13 Polyclonal Antibody

KPAB0346

Description

This Biotinylated Anti-human IL-13 Polyclonal Antibody has been validated for ELISA, ELISpot applications.

Product Information

SKU: KPAB0346
Conjugated: Biotin
Applications: ELISA, ELISpot

Antibody Data

Host: Rabbit
Reactivity: Human
Immunogen: Recombinant Human IL-13
Purification: Antigen-affinity Purification
Homology: Homo sapiens (human)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.