

## Anti-human IL-17A Polyclonal Antibody

**KPAB0347**

### Description

---

This Anti-human IL-17A Polyclonal Antibody has been validated for ELISA, ELISpot, FC, Neutralization, WB applications.

### Product Information

---

<b>SKU:</b>	KPAB0347
<b>Concentration:</b>	1.0 mg/mL
<b>Conjugated:</b>	Unconjugated
<b>Applications:</b>	ELISA, ELISpot, FC, Neutralization, WB

### Antibody Data

---

<b>Host:</b>	Rabbit
<b>Reactivity:</b>	Human
<b>Immunogen:</b>	Recombinant Human IL-17A
<b>Molecular Weight:</b>	15.1 kDa (calculated) for human IL-17A
<b>Purification:</b>	Antigen-affinity Purification
<b>Formulation:</b>	1X PBS buffer containing 0.09% NaN <sub>3</sub>
<b>Cross-Reactivity:</b>	Bovine IL-17A: Weak; Canine IL-17A: Strong; Cynomolgus IL-17A: Moderate; Dolphin IL-17A: Moderate; Equine IL-17A: Moderate; Ferret IL-17A: Moderate; Ovine IL-17A: Weak; Rabbit IL-17A: Weak
<b>Homology:</b>	Homo sapiens (human)

## Preparation & Storage

**Storage:** Stable for up to twelve months from date of receipt at 2-8°C.

**Shipping:** Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.