

Anti-human GM-CSF Polyclonal Antibody

KPAB0349

Description

This Anti-human GM-CSF Polyclonal Antibody has been validated for ELISA, ELISpot, FC, Neutralization, WB applications.

Product Information

SKU:	KPAB0349
Concentration:	1.0 mg/mL
Conjugated:	Unconjugated
Applications:	ELISA, ELISpot, FC, Neutralization, WB

Antibody Data

Host:	Rabbit
Reactivity:	Human
Immunogen:	Recombinant Human GM-CSF
Molecular Weight:	14.5 kDa (calculated) for human GM-CSF
Purification:	Antigen-affinity Purification
Formulation:	1X PBS buffer containing 0.09% NaN3
Cross-Reactivity:	Cynomolgus GM-Csf: Weak; Rabbit GM-Csf: Weak
Homology:	Homo sapiens (human)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.