

Biotinylated Anti-human GM-CSF Polyclonal Antibody

KPAB0350

Description

This Biotinylated Anti-human GM-CSF Polyclonal Antibody has been validated for ELISA applications.

Product Information

SKU: KPAB0350

Conjugated: Biotin

Applications: ELISA

Antibody Data

Host: Rabbit

Reactivity: Human

Immunogen: Recombinant Human GM-CSF

Purification: Antigen-affinity Purification

Homology: Homo sapiens (human)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.