

Anti-human LIF Polyclonal Antibody

KPAB0351

Description

This Anti-human LIF Polyclonal Antibody has been validated for ELISA, ELISpot, FC, Neutralization, WB applications.

Product Information

SKU:	KPAB0351
Concentration:	1.0 mg/mL
Conjugated:	Unconjugated
Applications:	ELISA, ELISpot, FC, Neutralization, WB

Antibody Data

Host:	Rabbit
Reactivity:	Human
Immunogen:	Recombinant Human LIF
Molecular Weight:	19.7 kDa (calculated) for human LIF
Purification:	Antigen-affinity Purification
Formulation:	1X PBS buffer containing 0.09% NaN ₃
Homology:	Gorilla gorilla gorilla (western lowland gorilla), Homo sapiens (human), Hylobates moloch (silvery gibbon), Nomascus leucogenys (northern white-cheeked gibbon), Pan paniscus (pygmy chimpanzee), Pan troglodytes (chimpanzee), Pongo abelii (Sumatran orangutan)
Cross-Reactivity:	Bovine LIF: Strong; Canine LIF: Strong; Caprine LIF: Strong; Cynomolgus LIF: Strong; Dolphin LIF: Strong; Feline LIF: Moderate; Guinea Pig LIF: Moderate; Mouse LIF: Weak; Rabbit LIF: Moderate; Rat LIF: Weak; Swine LIF: Moderate

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.