

Biotinylated Anti-human LIF Polyclonal Antibody

KPAB0352

Description

This Biotinylated Anti-human LIF Polyclonal Antibody has been validated for ELISA, ELISpot applications.

Product Information

SKU: KPAB0352

Conjugated: Biotin

Applications: ELISA, ELISpot

Antibody Data

Host: Rabbit

Reactivity: Human

Immunogen: Recombinant Human LIF

Purification: Antigen-affinity Purification

Homology: Gorilla gorilla gorilla (western lowland gorilla), Homo sapiens (human), Hylobates moloch (silvery gibbon), Nomascus leucogenys (northern white-cheeked gibbon), Pan paniscus (pygmy chimpanzee), Pan troglodytes (chimpanzee), Pongo abelii (Sumatran orangutan)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.