

Anti-cynomolgus monkey CXCL9 Polyclonal Antibody

KPAB0357

Description

This Anti-cynomolgus monkey CXCL9 Polyclonal Antibody has been validated for ELISA, ELISpot, FC, Neutralization, WB applications.

Product Information

SKU:	KPAB0357
Concentration:	1.0 mg/mL
Conjugated:	Unconjugated
Applications:	ELISA, ELISpot, FC, Neutralization, WB

Antibody Data

Host:	Rabbit
Reactivity:	Cynomolgus Monkey
Immunogen:	Recombinant Cynomolgus Monkey CXCL9
Molecular Weight:	11.8 kDa (calculated) for cynomolgus monkey CXCL9
Purification:	Antigen-affinity Purification
Formulation:	1X PBS buffer containing 0.09% NaN ₃
Cross-Reactivity:	Human CXCL9: Moderate
Homology:	Chlorocebus sabaeus (green monkey), Macaca fascicularis (crab-eating macaque), Macaca nemestrina (pig-tailed macaque), Mandrillus leucophaeus (drill), Theropithecus gelada (gelada)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.